

Skin Code^{PRO}

The Newest Generation- In
Advanced Microneedling



One device. All skin types.

All skin Resurfacing.

Face			
Fine Lines	Deep Wrinkles	Acne Scars	Tightening
			
Skin Reconstruction		Appearance of Large Pores	
			
Body			
Stretch Marks	Scars from injury or surgery		Loose & Sagging Skin
			





Collagen induction therapy, also known as skin needling, stands out as one of the safest skin resurfacing treatments available. During the procedure, a single-use needle gently penetrates the skin, stimulating the release of growth factors. These factors aid in the natural production of collagen and elastin in the treated area. Notably, this method of collagen induction preserves the epidermis and melanocytes in the deep layers of the skin, preventing the possibility of post-inflammatory hyperpigmentation. Unlike procedures such as lasers or dermabrasion, RF-powered microneedling is suitable for most skin types and colors, including sensitive areas near the eyes. It's a targeted treatment with minimal impact on surrounding skin and boasts a rapid recovery time.

Advance Vacuum

The updated RF microneedling system integrates vacuum technology to enhance topical penetration. With a unique dual air chamber design, it effectively pushes topicals deeper into the skin. Enhanced tissue engagement ensures optimal air capture, while microneedles create channels for precise delivery of RF energy.


As the needles retract, air is released toward the skin, facilitating deeper penetration of topicals into the dermis. This results in an impressive increase in topical penetration, reaching up to 67% with both M and F Type needles. This advancement takes RF microneedling to unprecedented levels, offering unparalleled versatility. The system allows for customized treatments tailored to individual patient needs, catering to all skin types and areas of the body, regardless of the season.

Auto Capture Skin

Crear channels


Release vacuum &
needle

Finished




Skin Code^{PRO}

2MHz RF Frequency



Metal Frame



Stepping Motor Driving



Button Switch



Modular Component




360-degree silent wheel




MICRO NEEDLE RF Applicator

- MICRO NEEDLE (RF) applicator for procedures requiring ablation and resurfacing of the skin. Add vacuum - suction combined probe for better skin contact, especially around eyes and neck areas.




Stepping Motor Type

Stepping Motor Type, different to existing solenoid type, the needle inserts to the skin smoothly with no shock, and causes no bleeding and no pain after the procedure



Gold Plating


Needle is durable and also have high biocompatibility by applying Gold Plating. Patient with Metal allergy could also use it with not concerning Contact Dermatitis.



Needle Depth Control

0.4-4.0mm [0.2mm Step]

Operate epidermis and dermis layer by controlling the needle depth in unit of 0.2mm.



Safety Needle System

Sterilized disposable needle tip operator can easily notice the applying RF energy from LED light.


Needle Thickness: 0.2mm

The needle structure is easy to penetrate the skin with minimum resistance.





INSULATED & NON INSULATED

Insulated




Non-insulated



FRACTIONAL RF Applicator


- Non-invasive wrinkle reduction
- Suitable for the patients with fear of needles
- No down time
- Quick Recovery



■ Eyes
■ Nose Area



Face Area




Body Area




MONO/BIPOLAR

Under Mono



Under Bipolar



■ NEGATIVE POLE


■ POSITIVE ELECTRODE

Negative Handle: Together with MONO polar, energy transfer Deeper

Cooling handle: Max Reduce Pain



Before & After





Platinum SkinTech Limited
info@platinumskintechco.uk

